



## **BA 91** **Dissection table**

(picture of cod. 10886, BA91 fixed dissection table)

**The BA91 is the only one available dissection table.**

### **Code 10886 – BA91 FIXED DISSECTION TABLE**

#### **MAIN FEATURES**

- All made in stainless steel AISI 304 with soldering without porosity
- Divided into two tanks sections, each one of cm. 50x40x20H where one is perforated to facilitate liquid drainage and formalin steam aspiration,
- Shaped lower zone for liquids gathering, inclined towards the exhausts

**SIZES & WEIGHT NET:** 1500Wx820Dx870/1270H mm, about 180 kg

**APPROX shipping sizes & weights:** 1700Wx1100Dx1500H mm, about 210 Kg

\*\*\*\*\*



## **EG series** **Autopsy tables**

(picture of cod. 10834, EG102 - raisable autopsy table)

**The EG table are all made in stainless steel AISI18/10 and are available with different equipment: see ahead for detailed specs per model**

## **EG TABLES: AVAILABLE MODELS AND ORDERING INFORMATION**

### **Code 10833 - EG 101 RAISING AND ROTATING AUTOPSY TABLE**

#### **STRUCTURE**

- Totally made with 18/10 AISI 304 stainless steel with continuous wire welding without porosity.
- Work table top inclined through the dumping.
- 40 mm. grid dumping.
- Anti-overflow perimetrical edge.
- Rounded corners for a better cleaning and disinfection.
- Tank head: mm. 500 x 400 x 200 H approx. equipped with dumping.
- Column base: 750 mm. fixable on the floor by four automatic kingpins.

#### **Provided Equipment:**

- Work table top with four shelves on the lower washing level.
- Hot and cold water with mobile drain nozzle fixed on central column
- Mobile drainpipe near the tank head.
- Washing system with perimetrical nozzles and manual system command.
- Liquids aspiration system by Venturi's principle water powered (pressure higher than 2 bar).
- Smoke-aspiration and air-change managed by a centrifugal motor. Air-flow: 2000 mc/h. command switch is on the control panel.
- The centrifugal motor is equipped with an anchorage stirrup for external application.
- Electric monophasic plug: 10 A.
- Adjustable top (max: 180 mm.)

#### **MOVEMENTS**

- 90° rotation on both left and right side from the central position: total rotation 180°
- movement by hydraulic plant integrated in the central column (low noise level, electric motor with 24-volt voltage continuous current)
- movements controls by a four-direction master controller frontal- positioned for easy use.

**SIZES & WEIGHT NET:** 2720W x 900D x 850H mm, about 250 kg

**APPROX shipping sizes & weights:** 2950Wx1150Dx1300H mm, about 270 kg

\*\*\*\*\*

### **Code 10834 - EG 102 RAISING ONLY AUTOPSY TABLE**

#### **STRUCTURE**

- Totally made with 18/10 AISI 304 stainless steel with continuous wire welding without porosity.
- Work table top inclined through the dumping.
- 40 mm. grid dumping.
- Anti-overflow perimetrical edge.
- Rounded corners for a better cleaning and disinfection.
- Tank head: mm. 500 x 400 x 200 H approx. equipped with dumping.

- Column base: 750 mm. fixable on the floor by four automatic kingpins.

### Provided Equipment:

- Work table top with four shelves on the lower washing level.
- Hot and cold water with mobile drain nozzle fixed on central column
- Mobile drainpipe near the tank head.
- Washing system with perimetrical nozzles and manual system command.
- Liquids aspiration system by Venturi's principle water powered (pressure higher than 2 bar).
- Smoke-aspiration and air-change managed by a centrifugal motor. Air-flow: 2000 mc/h. command switch is on the control panel.
- The centrifugal motor is equipped with an anchorage stirrup for external application.
- Electric monophasic plug: 10 A.
- Adjustable top (max: 180 mm.)

### MOVEMENTS

- Raisable range up – down up to 180 mm
- movement by hydraulic plant integrated in the central column (low noise level, electric motor with 24-volt voltage continuous current)

**SIZES & WEIGHT NET:** 2720W x 900D x 850H mm, about 230 kg

**APPROX shipping sizes & weights:** 2950Wx1150Dx1300H mm, about 260 kg

\*\*\*\*\*

## Code 10835 - EG 103 FIXED AUTOPSY TABLE

### STRUCTURE

- Totally made with 18/10 AISI 304 stainless steel with continuous wire welding without porosity.
- Work table top inclined through the dumping.
- 40 mm. grid dumping.
- Anti-overflow perimetrical edge.
- Rounded corners for a better cleaning and disinfection.
- Tank head: mm. 500 x 400 x 200 H approx. equipped with dumping.
- Column base: 750 mm. fixable on the floor by four automatic kingpins.

### Provided Equipment:

- Work table top with four shelves on the lower washing level.
- Hot and cold water with mobile drain nozzle fixed on central column
- Mobile drainpipe near the tank head.
- Washing system with perimetrical nozzles and manual system command.
- Liquids aspiration system by Venturi's principle water powered (pressure higher than 2 bar).
- Smoke-aspiration and air-change managed by a centrifugal motor. Air-flow: 2000 mc/h. command switch is on the control panel.
- The centrifugal motor is equipped with an anchorage stirrup for external application.
- Electric monophasic plug: 10 A.
- Adjustable top (max: 180 mm.)

## MOVEMENTS

- None

**SIZES & WEIGHT NET:** 2720W x 900D x 850H mm, about 220 kg

**APPROX shipping sizes & weights:** 2950Wx1150Dx1300H mm, about 255 kg

\*\*\*\*\*

## Code 10836 - EG 104 BASIC AUTOPSY TABLE

### STRUCTURE

- Totally made with 18/10 AISI 304 stainless steel with continuous wire welding without porosity.
- Work table top inclined through the dumping.
- 60 mm. grid dumping.
- Anti-overflow perimetrical edge.
- Rounded corners for a better cleaning and disinfection.
- Tank head: mm. 290x400x 200H mm approx. equipped with dumping.
- Base with tubular legs 30x30mm, fixable on the floor by four automatic kingpins.

### Provided Equipment:

- Work table top with four shelves on the lower washing level.
- Hot and cold water with mobile drain nozzle near the tank head
- Fixed top at 850 mm

## MOVEMENTS

- None

**SIZES & WEIGHT NET:** 2600W x 750D x 850H mm, about 215 Kg

**APPROX shipping sizes & weights:** 2800Wx950Dx1100H mm, about 240 kg

\*\*\*\*\*

## OPTIONAL ACCESSORIES

Code	Description
10563	Table for organs carrying for EG101, 102 and 103 only
11485	Table for organs carrying for EG104 only
11766	Charcoal filtering system for EG101, 102, 103 and for dissection table BA91
10872	Teflon headrest

## **OTHERS TECHNICAL DATA**

**Product Names :**

**EG autopsy tables  
BA91 dissection table**

**Manufacturer:**

Angelantoni Life Science S.r.L.  
Loc. Cimacolle, 464  
06056 Massa Martana (PG) - Italy

**Ambient conditions of use:**

Temperature: +10° ~ +32° C  
Relative humidity: 30 ~ 80% without condensation

**Stocking condition:**

Temperature: +2° ~ + 50° C  
Relative humidity: 20 ~ 80% without condensation

**Alarms types:**

not applicable

**In case of power black-out:**

all electrical equipment stop