



# Pass Box


Pass-Through Chamber  
Material Transfer Hatch

- Stainless steel interior & exterior housing with mirror polish finish
  - Large radius corners making it perfectly cleanable & sealable
  - Power supply free innovative reliable interlock (mechanical)
    - Installation parts included as standard
    - Variety of accessories and dimensions
      - Tempered 12mm glass door
      - More than 500 installations
      - Single-handed operation





### Characterization



### Most Frequent Dimensions / Models

	Interior			Exterior			Equipment	View
	W	H	D	W	H	D		
<b>PB 4/4/4</b>	400	400	400	500	500	450	/	
<b>PB 5/5/5</b>	500	500	500	600	600	550	/	
<b>PB 6/6/6</b>	600	600	600	700	700	650	/	
<b>PB 8/8/8</b>	800	800	800	900	900	850	/	
<b>PB 6/6/6 UV</b>	600	600	600	700	850	750	UV – C surface sterilization	
<b>PB 6/6/6 HUV</b>	600	600	600	750	970	750	Internal air circulation through HEPA filter sterilization	
<b>PB 6/6/6 AIR</b>	600	600	600	700	700	750	Air shower with supply from HVAC system	
<b>PB 6/6/6 EL</b>	600	600	600	700	750	750	Electrical interlock	
<b>PB 8/10/8 LUX</b>	800	1000	800	1000	1150	950	Floor installation	
<b>PB 6/6/6 90°</b>	600	600	600	*	*	*	Corner installation	

All other dimensions available on demand!

